



POLYPURE® S+

High Flux Series





POLYPURE® S+



Allmed's high flux dialyser, superbly balancing middle molecule removal, optimal albumin retention, whilst offering excellent levels of endotoxin protection consistent with POLYPURE® membranes. All feature the additional advantages of autoclaved steam sterilisation.

Steam sterilisation in autoclave: best for patients & the environment

- No potentially dangerous substances are generated during autoclaved steam sterilisation
- Lower priming volume requirements compared to other sterilisation methods for the same product (see DFU)
- Optimised rheology during clinical use
- Autoclave steam sterilisation is the only steam sterilisation recognised by The International Pharmacopeia (Methods of Sterilisation-VI Edition)
- Less impact on the environment compared with irradiation and ethylene oxide

Excellent balance between middle molecule removal and albumin retention: safeguarding the patient nutritional status

The membrane structure perfects the delicate balance between substantial middle molecule removal whilst preserving levels of albumin. POLYPURE® S+ removes 85% of β_2 microglobulin with only 0.6g of albumin loss, during a session of post dilution HDF.



Micro-undulated technology: highest efficacy

Micro-undulation technology ensures the highest efficiency and clearances thanks to “bundles” and “dead-zones” prevention, together with the optimally increased surface exchange of blood-membrane-dialysate.


Protection from endotoxins: reduced inflammation

POLYPURE® asymmetrical membrane with an optimised micro-undulated fibre structure and 40µm thickness, enhances the membrane’s ability to prevent endotoxins transfer due to its geometrical and adsorptive capacity capability.

Systemic inflammation is well recognised to be associated with an increased cardiovascular disease risk in patients with chronic renal disease. Exposure to endotoxins results in a release of a wide variety of pro-inflammatory cytokines, with subsequent chronic systemic inflammation and a broad range of negative cardiovascular effects. Chronic inflammation can also modify the process of erythropoiesis and probably induces resistance to erythropoietin.



Allmed Laboratories–Electron Microscopy
Image of Allmed membrane



IGMHS 20080220rdc x 5000 10µm

POLYPURE® S+		10 S+	13 S+	14 S+	16 S+	18 S+	20 S+	22 S+
Q _B = 200 ml/min	Urea	185	191	192	193	194	196	197
	Creatinine	170	178	180	183	187	189	191
	Phosphate	158	172	176	178	183	186	188
	Vitamin B12	108	126	133	136	145	148	153
	Inulin	74	89	93	97	109	115	121
Q _B = 300 ml/min	Urea	241	258	261	264	270	275	279
	Creatinine	207	225	231	234	245	250	256
	Phosphate	191	211	219	225	236	241	247
	Vitamin B12	121	144	152	157	170	175	183
	Inulin	80	98	103	110	121	129	137
Q _B = 400ml/min	Urea	275	299	303	310	319	328	335
	Creatinine	231	254	261	267	281	289	298
	Phosphate	211	237	243	254	270	278	287
	Vitamin B12	129	155	163	170	186	191	201
	Inulin	84	103	108	115	129	136	145
Clearance in vitro (ml/min)								
QD = 500 ml/min		Surface area (m2)	1.0	1.3	1.4	1.6	1.8	2.2
QF = 0 ml/min		UF coefficient (ml/hr*mmHg)	32	43	47	53	58	74
T = 37°C		Blood priming Volume (ml)	59	69	75	86	105	115
Max TMP = 500mmHg			▼	▼	▼	▼	▼	▼
Sieving coefficient	Vitamin B12	1.0						
	Inulin	1.0						
	β2 microglobulin	0.85						
	Albumin	< 0.001						

- Performance data was measured in vitro according to standard ISO 8637
- UF measurement: using bovine/human blood (Hct 32%;protein 60g/l)
- Typical values obtained with an individual batch of fibres, clinical use may illustrate a difference in results in relation to different Ultrafiltration/measuring techniques and possible variation between batches of fibres

DISCLAIMER: This document was produced November 2025 to provide a general outline only. Information, illustrations, images or figures (together "materials") contained in this document (which does not purport to be comprehensive) are illustrative and/or indicative only and are subject to change.

While every effort has been made to ensure the accuracy of this publication, [Allmed Medical GmbH], KIMAL plc and its group of companies accept no liability in relation to the adequacy, accuracy, completeness or reasonableness of the materials. All and any such responsibility and liability are expressly disclaimed to the fullest extent permitted by law.



Contact us for more information:

Tel: **+44 (0)845 437 9542**

Email: renal@kimal.com

Scan the QR code
to learn more



www.kimal.com



BRO-GEN-(18) Kimal | Allmed - POLYPURE S+ | November 2025