

➤ KIMAL GUIDEWIRES: GREATER CONTROL FOR SELDINGER INSERTION



- Grooved thumb wheel for controlled advancement
- Available in sizes from .018" to .035"
- Kink-resistant J-tip nitinol wire
- Clear depth markings on both guidewire and dispenser

KIMAL HAS A SIMPLE SOLUTION TO DROPPED OR DAMAGED GUIDEWIRES...


REPLACEMENT GUIDEWIRES

All Kimal guidewires are nitinol to minimise kinking and are preloaded in the thumb wheel dispenser to help obtain vascular access.

The tip of the guidewire dispenser is designed to fit into the Altius® guidewire syringe for CVC placement.

Grooved thumb wheel provides greater control when advancing whilst wearing gloves. The anti-slip groove ensures the guidewire stays centred under your thumb.

Depth markings are indicated both on the guidewire dispenser and the guidewire.



REPLACEMENT GUIDEWIRES (SUPPLIED STERILE) PRODUCT ORDERING DETAILS

PRODUCT CODE	BOX QTY	MATERIAL	TIP	GAUGE	LENGTH
KGW01850N	25	Nitinol	J-tip	0.018" (0.46mm)	50cm
KGW02150N	25	Nitinol	J-tip	0.021" (0.53mm)	50cm
KGW03550N	25	Nitinol	J-tip	0.035" (0.89mm)	50cm
KGW03570N	25	Nitinol	J-tip	0.035" (0.89mm)	70cm

For more information, pricing or samples, speak to your local Kimal Representative, or contact:

UK
Tel: +44 (0)845 437 9542
Email: customer@kimal.com

INTERNATIONAL
Tel: +44 (0)845 437 9543
Email: internationalenquiries@kimal.com

MIDDLE EAST, ASIA & AFRICA
Tel: +97 16 552 8146
Email: mea@kimal.com