


SPECIALITY PRODUCTS

Radial Compression Band

Precision compression where it matters most



A radial compression band is intended for use following trans radial artery access procedures to support controlled haemostasis at the wrist puncture site. The device is designed to apply adjustable, localized pressure to assist in achieving haemostasis while supporting arterial patency.

The radial compression band is used as part of standard post-procedural care to help manage bleeding at the access site and to support the reduction of procedure-related complications.



DESIGNED FOR CONTROLLED COMPRESSION

The radial compression band is designed to exert consistent, adjustable pressure at the access site to support haemostasis and assist in managing localised swelling following trans radial procedures. The device enables clinicians to tailor compression according to individual patient needs and clinical protocols.

WHERE PRECISION MEETS PATIENT COMFORT

- + Support for Patent Haemostasis:** The device is designed to allow the application of controlled pressure intended to support haemostasis while maintaining arterial blood flow, in line with established clinical practice.
 - + Support for Complication Management:** By providing localised, adjustable compression at the puncture site, the radial compression band is intended to support post-procedural care and help reduce the risk of access-site complications, such as hematoma formation, when used appropriately.
 - + Patient Comfort:** An inflatable, transparent balloon enables gradual adjustment of compression, supporting patient comfort during use.
 - + Monitoring Capability:** The transparent design allows visual observation of the access site for signs of bleeding or bruising while the device is in place.
 - + Fit and Adjustability:** Two band lengths are available to support appropriate fit across a range of patient anatomies.
-

Our **difference**
is your **speciality**

ORDERING INFORMATION

Product Code	Band Type	Length (cm)
42.17.30001	Regular	24
42.17.30004	Large	27

Refer to the information for use and package insert for precautions, warnings, potential complications, and full instructions for use.

For more information, pricing or samples speak to your local Kimal Clinical Specialist or contact:

UK

Tel: +44 (0)845 437 9542
Fax: +44 (0)845 437 9544
Email: customer care@kimal.com

International

Tel: +44 (0)845 437 9543
Fax: +44 (0)845 437 9549
Email: internationalenquiries@kimal.com

France

Tel: +33 (0)474 13 74 23
Fax: +33 (0)4 74 13 70 13
Email: info@kimal.fr

© 2026 Kimal PLC. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal PLC.

Kimal PLC reserve the right to change the specification of its products at any time.
Kimal PLC, Clayfield Road, Worcester Six Business Park, Worcestershire WR4 0AE.